

## Product Datasheet

# Human Neuregulin-1 protein (Pro-Ser2-Lys246) (orb311953)

<b>Description</b>	Recombinant Human Neuregulin-1 beta is produced by our E.coli expression system and the target gene ...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	tag free
<b>Note</b>	For research use only
<b>Uniprot ID</b>	AAA58639.1
<b>Expiration Date</b>	6 months from date of receipt.