

## Product Datasheet

### Human PACSIN1 protein (C-His) (orb311963)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Recombinant Human Protein Kinase C and Casein Kinase Substrate in Neurons Protein 1 is produced by o... |
| <b>Endotoxins</b>        | Less than 0.1 ng/ug (1 EU/ug)   |
| <b>Form/Appearance</b>   | Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.                                  |
| <b>Storage</b>           | Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.                  |
| <b>Tag</b>               | C-His   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | 6His tag at the C-terminus.   |
| <b>Uniprot ID</b>        | <a href="#">Q9BY11</a>  |
| <b>Expiration Date</b>   | 6 months from date of receipt.  |