

Product Datasheet

Rat MCSF protein (orb311970)

Description	Recombinant Rat Macrophage colony-stimulating factor 1 is produced by our Mammalian expression syste...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	Q8JZQ0
Expiration Date	6 months from date of receipt.