

## Product Datasheet

# Human Receptor-Type Tyrosine-Protein Phosphatase-Like N protein (N-His)(Arg576-Gln950) (orb311972)

<b>Description</b>	Recombinant Human Islet Cell Antigen 2/Protein Tyrosine Phosphatase Receptor-type N is produced by o...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	N-His
<b>Note</b>	For research use only
<b>Application notes</b>	6His tag at the N-terminus.
<b>Uniprot ID</b>	<a href="#">Q16849</a>
<b>Expiration Date</b>	6 months from date of receipt.