

Product Datasheet

Human Receptor-Type Tyrosine-Protein Phosphatase-Like N protein (N-His)(Met687-Gln979) (orb311974)

Description	Recombinant Human Islet Cell Antigen 2/Protein Tyrosine Phosphatase Receptor-type N is produced by o...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Uniprot ID	Q16849
Expiration Date	6 months from date of receipt.