

Product Datasheet

Human S100A8 protein (C-His) (orb311986)

Description	Recombinant Human S100 Calcium Binding Protein A8 is produced by our E.coli expression system and th...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mMHEPES,150mMnacl,2.5mMEDTA,10% Glycerol,ph7.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	P05109
Expiration Date	6 months from date of receipt.