

Product Datasheet

Human MECR protein (C-His) (orb312045)

Description	Recombinant Human Trans-2-Enoyl-CoA Reductase Mitochondrial is produced by our Mammalian expression ...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	AAH01419.1
Expiration Date	6 months from date of receipt.