

Product Datasheet

Human Transcriptional Repressor Protein YY1 protein (C-His) (orb312047)

Description	Recombinant Human Yin and Yang 1 protein/Transcriptional repressor protein YY1 is produced by our E....
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	P25490
Expiration Date	6 months from date of receipt.