

Product Datasheet

Human TGF Beta 3 protein (orb312049)

Description	Recombinant Human/Mouse Transforming Growth Factor beta 3 is produced by our Mammalian expression sy...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 50mM Glycine-HCl, 150mM NaCl, pH 2.5.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	tag free
Note	For research use only
Uniprot ID	P10600
Expiration Date	6 months from date of receipt.