

Product Datasheet

Mouse CD262 protein (C-Fc) (orb312057)

Description	Recombinant Mouse TNF-Related Apoptosis-Inducing Ligand Receptor 2 is produced by our Mammalian expr...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc
Note	For research use only
Application notes	Fc tag at the C-terminus.
Uniprot ID	Q9QZM4
Expiration Date	6 months from date of receipt.