

Product Datasheet

Human PTP4A2 protein (C-His) (orb312068)

Description	Recombinant Human Phosphatase of Regenerating Liver 2 is produced by our E.coli expression system an...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of 20mM Hepes,150mMNacl,10mMβ-ME,pH7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	Q12974
Expiration Date	6 months from date of receipt.