

Product Datasheet

phosphor-c-Abl (phospho-Tyr412) antibody (orb312208)

Description	phosphor-c-Abl (phospho-Tyr412) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human c-Abl (LA(p-Y)NK)
Target	phosphor-c-Abl(Tyr412)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	P00519
Expiration Date	12 months from date of receipt.