

Product Datasheet

FMDV VP1 antibody (orb312248)

Description	FMDV VP1 antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant FMDV Asia1 VP1 (His tag)
Target	FMDV VP1
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	240 kDa
Uniprot ID	P03305
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.