

Product Datasheet

MAP2 (phospho-Thr1613/1619) antibody (orb312314)

Description

MAP2 (phospho-Thr1613/1619) antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

KLH conjugated synthesised phosphopeptide derived from human MAP2 around the phosphorylation site of Thr1613 and Thr1619 (SR(p-T)PGTPG(p-T)PS)

Target

Phospho-MAP2 (Thr1613/1619)

Preservatives

Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

70/201 kDa

Uniprot ID


[P11137](#)

Dilution Range


WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=2ug/Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.



Flow
cytometry of
SH-
SY5Y.Primary
Antibo...



Western
blot:Lane
1:SH-SY5Y
(Human)
Cell...