

Product Datasheet

NFkB p100 (phospho-Ser869) antibody (orb312388)

Description	NFkB p100 (phospho-Ser869) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NFkB p100 around the phosphorylation site of Ser870 (YG(p-S)QS)
Target	phospho-NFkB2 (Ser870)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50/100 kDa
Uniprot ID	Q00653
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.