

Product Datasheet

SSX4 antibody (orb31269)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target SSX4

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 21858

Uniprot ID [O60224](#)

NCBI [NP_001035702.1](#), [NP_783856.1](#), [NP_001030004.1](#), [NP_005627.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

