

Product Datasheet

PKC zeta (phospho-Thr560) antibody (orb312763)

Description	PKC zeta (phospho-Thr560) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PKC zeta around the phosphorylation site of Thr560 (QL(p-T)PD)
Target	phospho-PKC zeta (Thr560)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	67 kDa
Uniprot ID	Q05513
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.