

Product Datasheet

PLB(phospho-Ser16+Thr17) antibody (orb312832)

Description	PLB(phospho-Ser16+Thr17) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PLB around the phosphorylation site of Ser16+Thr17 (RA(p-S)(p-T)IE)
Target	phospho-PLB (Ser16+Thr17)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6 kDa
Uniprot ID	P26678
Expiration Date	12 months from date of receipt.