

Product Datasheet

Survivin (phospho-Thr117) antibody (orb313219)

Description	Survivin (phospho-Thr117) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Survivin around the phosphorylation site of Thr117 (KE(p-T)NN)
Target	phospho-Survivin (Thr117)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16.5 kDa
Uniprot ID	O15392
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.