

Product Datasheet

TK1 (phospho-Ser13) antibody (orb313258)

Description	TK1 (phospho-Ser13) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Thymidine Kinase 1 around the phosphorylation site of Ser13 (PG(pS)PS)
Target	phospho-Thymidine Kinase 1 (Ser13)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	P04183
Expiration Date	12 months from date of receipt.