

Product Datasheet

SHB (phospho-Tyr246) antibody (orb313260)

Description	SHB (phospho-Tyr246) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human SHB around the phosphorylation site of Tyr246 (Y(p-S)DP)
Target	phospho-SHB (Tyr246)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	55 kDa
Uniprot ID	Q15464
Expiration Date	12 months from date of receipt.