


Product Datasheet

Human ADAMTS4 ELISA Kit (orb314901)

Description	Sandwich ELISA for quantitative detection of Human ADAMTS4 concentrations.
Reactivity	Human
Range	625 pg/ml - 40000 pg/ml
Tested Applications	ELISA
Target	ADAMTS4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	075173
Sensitivity	20 pg/ml
Entrez	9507
Expiration Date	Please enquire.



Representative
Standard
Curve