


Product Datasheet

Human Surfactant Protein D ELISA Kit (orb314929)

Description	Sandwich ELISA for quantitative detection of Human Surfactant Protein D concentrations.
Reactivity	Human
Range	625 pg/ml - 40000 pg/ml
Tested Applications	ELISA
Target	SFTPD
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P35247
Sensitivity	20 pg/ml
Entrez	6441
Expiration Date	Please enquire.



Representative Standard Curve