


## Product Datasheet

### Mouse GPNMB ELISA Kit (orb314950)

<b>Description</b>	Sandwich ELISA for quantitative detection of Mouse GPNMB concentrations.
<b>Reactivity</b>	Mouse
<b>Range</b>	156 pg/ml - 10000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	GPNMB
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma
<b>Uniprot ID</b>	<a href="#">Q99P91</a>
<b>Sensitivity</b>	10 pg/ml
<b>Entrez</b>	<a href="#">93695</a>
<b>Expiration Date</b>	Please enquire.



Representative  
Standard  
Curve