


## Product Datasheet

### Mouse TIG2 ELISA Kit (orb314962)

<b>Description</b>	Sandwich ELISA for quantitative detection of Mouse TIG2 concentrations.
<b>Reactivity</b>	Mouse
<b>Range</b>	780 pg/ml - 50000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	RARRES2
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma
<b>Uniprot ID</b>	<b>Q9DD06</b>
<b>Sensitivity</b>	20 pg/ml
<b>Entrez</b>	<b>71660</b>
<b>Expiration Date</b>	Please enquire.



Representative  
Standard  
Curve