


Product Datasheet

Mouse Calprotectin ELISA Kit (orb314963)

Description	Mouse Calprotectin ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	S100A9
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P31725
Sensitivity	10 pg/ml
Entrez	20202
Expiration Date	Please enquire.



Representative Standard Curve