


Product Datasheet

Human Resistin ELISA Kit (orb315009)

Description	Sandwich ELISA for quantitative detection of Human Resistin concentrations.
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	RETN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9HD89
Sensitivity	4 pg/ml
Entrez	56729
Expiration Date	Please enquire.



Representative
Standard
Curve