


Product Datasheet

Human IL-17 ELISA Kit (orb315014)

Description	Sandwich ELISA for quantitative detection of Human IL-17 concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	IL17A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q16552
Sensitivity	7 pg/ml
Entrez	3605
Expiration Date	Please enquire.



Representative Standard Curve