


Product Datasheet

Human Visfatin ELISA Kit (orb315023)

Description	Sandwich ELISA for quantitative detection of Human Visfatin concentrations.
Reactivity	Human
Range	1000 pg/ml - 64000 pg/ml
Tested Applications	ELISA
Target	NAMPT
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P43490
Sensitivity	500 pg/ml
Entrez	10135
Expiration Date	Please enquire.



Representative
Standard
Curve