


Product Datasheet

Mouse IL-17 ELISA Kit (orb315053)

Description	Sandwich ELISA for quantitative detection of Mouse IL-17 concentrations.
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL17A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	Q62386
Sensitivity	15 pg/ml
Entrez	16171
Expiration Date	Please enquire.



Representative Standard Curve