


## Product Datasheet

### Pig IL-2 ELISA Kit (orb315066)

<b>Description</b>	Sandwich ELISA for quantitative detection of Pig IL-2 concentrations.
<b>Reactivity</b>	Porcine
<b>Range</b>	3.9 pg/ml - 250 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL2
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
<b>Sensitivity</b>	2 pg/ml
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve