


Product Datasheet

Mouse TNFsR I ELISA Kit (orb315075)

Description	Mouse TNFsR I ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L30-A212
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum is 460-1890 pg/ml
Sample Types	cell culture supernatants and serum.
MW	50130 MW
Uniprot ID	P25118
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
TNFR1
PicoKine
ELISA Kit
standard ...