


Product Datasheet

Mouse CXCL9/Mig ELISA Kit (orb315082)

Description	Mouse CXCL9/Mig ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: T22-T126
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 0-710g/ml
Sample Types	cell culture supernatants, serum and plasma (EDTA).
MW	14444 MW
Uniprot ID	P18340
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
CXCL9
PicoKine
ELISA Kit
standard ...