


Product Datasheet

Human ADAMTS13 ELISA Kit (orb315094)

Description	Human ADAMTS13 ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: Q34-W688
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is around 500 ng/ml
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	153604 MW
Uniprot ID	Q76LX8
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human
ADAMTS13
PicoKine
ELISA Kit
standa...