


Product Datasheet

Mouse S100A8/Calgranulin A ELISA Kit (orb315096)

Description	Mouse S100A8/Calgranulin A ELISA Kit
Reactivity	Mouse
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: P2-E89
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10295 MW
Uniprot ID	P27005
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Mouse S100A8 PicoKine ELISA Kit standard...