


Product Datasheet

Human Relaxin 1 ELISA Kit (orb315105)

Description	Human Relaxin 1 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R162-Cys185 (A) & K26-S53 (B)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 78 pg/ml
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	21146 MW
Uniprot ID	P04808
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Relaxin 1 PicoKine ELISA Kit stand...