


Product Datasheet

Human Galectin-3BP ELISA Kit (orb315106)

Description	Human Galectin-3BP ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: V19-D585
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:2000, concentration in serum and plasma is around 1.04-22.08 ug/ml
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	65331 MW
Uniprot ID	Q08380
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Galectin-3BP PicoKine ELISA Kit st...