


Product Datasheet

Human Complement Factor H/CFH ELISA Kit (orb315117)

Description	Human Complement Factor H/CFH ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Target protein isolated from human plasma
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:20000, concentration in serum and plasma is around 300 ug/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	139096 MW
Uniprot ID	P08603
Sensitivity	50 pg/ml
Expiration Date	Please enquire.



Human Complement H/CFH PicoKine ELISA Ki...