

Product Datasheet

Hsp90 alpha/HSP90AA1 Antibody (orb315143)

Description

Hsp90 alpha/HSP90AA1 Antibody

Species/Host

Rabbit

Reactivity

Human, Monkey, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human Hsp90 alpha (454-488aa QNRKKLSELLRYYTSASGDEMVSCLKDYCTRMKEN Q), identical to the related mouse and rat sequences.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


84660 MW

Uniprot ID


P07900

Expiration Date


12 months from date of receipt.



Flow Cytometry analysis of SiHa cells us...



WB analysis of Hsp90 alpha using anti-Hs...



IF analysis of Hsp90 alpha using anti-Hs...