


## Product Datasheet

### SA4503 2HCl (orb315320)

<b>Description</b>	SA4503 is a selective sigma 1 receptor( $\sigma$ 1R) agonist; high affinity for the sigma 1 receptor subtype...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{23}H_{34}Cl_2N_2O_2$
<b>MW</b>	441.437
<b>CAS Number</b>	[165377-44-6]
<b>SMILES</b>	<chem>COC1=C(C=C(C=C1)CCN2CCN(CC2)CCCC3=CC=CC=C3)OC.Cl.Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.



Chemical structure of SA4503 2HCl