

## Product Datasheet

### TP-808 (orb315452)

<b>Description</b>	TP808 is a useful intermediate for the synthesis of diverse tetracycline antibiotics.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{26}H_{34}N_2O_5Si$
<b>MW</b>	482.64
<b>CAS Number</b>	[852821-06-8]
<b>SMILES</b>	<chem>CC(C)(C)[Si](C)(C)OC12C(CC=CC1=O)C(C3=C(C2=O)C(=NO3)OCC4=C C=CC=C4)N(C)C</chem>
<b>Expiration Date</b>	12 months from date of receipt.