


Product Datasheet

CD8a protein (orb315535)

Description	CD8a protein
Species/Host	E. coli
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, SDS-PAGE, WB
Target	CD8a
Form/Appearance	20mM Tris-HCl, 50mM NaCl, 1M Urea, 1mM DTT pH 7.5
Concentration	1 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Application notes	N-His6
Protein Sequence	SQFRVSPLDR TWNLGETVEL KCQVLLSNPT SGCSWLFQPR GAAASPTFLL YLSQNKPKAA EGLDTQRFSG KRLGDTFVLT LSDFRRENEG YYFCSALSNS IMYFSHFVPV FLPKPTTTP A
Purity	> 85.0% as determined by 15% SDS-PAGE
MW	17.0 kDa
NCBI	001768
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of CD8a protein