

Product Datasheet

Kir3.1 (Phospho-S185) antibody (orb315577)

Description Rabbit polyclonal antibody to Kir3.1 (Phospho-S185)

Species/Host Rabbit

Reactivity Human, Mouse, Primate, Rat

Conjugation Unconjugated

Tested Applications IF, IH, WB

Immunogen KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S185 of human Kir3.1 protein. The exact sequence is proprietary.

Target KCNJ3

Preservatives Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal


Source Rabbit

Uniprot ID [P48549](#), [P63250](#), [P63251](#)


Entrez [3760](#), [16519](#), [50599](#)

Dilution Range WB: 1-500:1000, IHC-P: 1-100:200, IF/ICC: 1-100:500


Expiration Date 12 months from date of receipt.



Western blot analysis of Kir3.1 (Phospho...



Immunohistochemical analysis of Kir3.1 (...)



Immunofluorescent analysis of Kir3.1 (Ph...