

Product Datasheet

c-Jun (Phospho-S63) peptide (orb315839)

Description	The peptide is used to c-Jun (pS63) antibody
Reactivity	Human, Mouse, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	JUN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05627 , P17325 , P05412
Entrez	16476 , 3725 , 24516
Expiration Date	12 months from date of receipt.