

Product Datasheet

c-Jun (Phospho-S73/100) peptide (orb315842)

Description	The peptide is used to c-Jun (pS73/100) antibody
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	JUN; JUND
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15066 , P05627 , P17325 , P05412 , P17535 , P52909
Entrez	16476 , 24518 , 16478 , 3725 , 3728 , 24516
Expiration Date	12 months from date of receipt.