

## Product Datasheet

# JNK1/2/3 (Phospho-T183/Y185) peptide (orb315867)

<b>Description</b>	The peptide is used to JNK1/2/3 (pT183/Y185) antibody
<b>Reactivity</b>	Human, Mouse, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	MAPK8; MAPK9; MAPK10
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q61831</a> , <a href="#">Q91Y86</a> , <a href="#">P49186</a> , <a href="#">P53779</a> , <a href="#">P49187</a> , <a href="#">P45983</a> , <a href="#">P45984</a> , <a href="#">Q9WTU6</a> , <a href="#">P49185</a>
<b>Entrez</b>	<a href="#">50658</a> , <a href="#">5602</a> , <a href="#">26419</a> , <a href="#">26420</a> , <a href="#">5599</a> , <a href="#">5601</a> , <a href="#">25272</a> , <a href="#">116554</a>
<b>Expiration Date</b>	12 months from date of receipt.