

Product Datasheet

Ezrin/Radixin/Moesin (Phospho-T567/564/558) peptide (orb315886)

Description	The peptide is used to Ezrin/Radixin/Moesin (pT567/564/558) antibody
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	MSN; RDX; EZR
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35241 , P15311 , P31977 , P26038 , O35763 , P26041 , P26043 , P26040
Entrez	17698 , 4478 , 54319 , 81521 , 22350 , 5962 , 19684
Expiration Date	12 months from date of receipt.