

Product Datasheet

PEA15 (Phospho-S104) peptide (orb315914)

Description	The peptide is used to PEA15 (pS104) antibody
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	PEA15
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q15121 , Q5U318 , Q62048
Entrez	364052 , 8682 , 18611
Expiration Date	12 months from date of receipt.