

Product Datasheet

MYPT1 (Phospho-T853) peptide (orb315923)

Description	The peptide is used to MYPT1 (pT853) antibody
Reactivity	Human, Mouse, Rat, Zebrafish
Tested Applications	Blocking
Target	PPP1R12A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O14974 , Q10728 , Q9DBR7
Entrez	17931 , 4659 , 116670
Expiration Date	12 months from date of receipt.