

Product Datasheet

PKA C alpha/beta/gamma (Phospho-T198) peptide (orb315925)

Description	The peptide is used to PKA C alpha/beta/gamma (pT198) antibody
Reactivity	Human, Mouse, Porcine, Rat, Sheep
Tested Applications	Blocking
Target	PRKACA; PRKACB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P68181 , P68182 , P22612 , P05132 , P22694 , P17612 , P27791
Entrez	18747 , 5567 , 5566 , 18749 , 293508
Expiration Date	12 months from date of receipt.